MouseDB Documentation

Release 0.1

Dave Bridges, Ph.D.

CONTENTS

1	MouseDB Concepts 1 Animal Module			3 3
2	.2 Data Module	• •	•	3 5
_	.1 Configuration			5
	.2 Software Dependencies			5
	.3 Database Setup			5
	.4 Web Server Setup			6
	.5 Final Configuration and User Setup		•	6
3	Animal Data Entry			7
	.1 Newborn Mice or Newly Weaned Mice			7
	.2 Newborn Mice			7
	.3 Weaning Mice			7
	.4 Cage Changes (Not Weaning)			7
	.5 Genotyping or Ear Tagging			8
	.6 Marking Mice as Dead		•	8
4	studies and Experimental Setup			9
5	Measurment Entry			11
	.1 Studies			11
	.2 Experiment Details			11
	.3 Measurements		•	11
6	Automated Documentation			13
	.1 Data Package			13
	.2 Animals Package			17
	.3 Timed Mating Package			20
	.4 Groups Package		•	20
7	ndices and tables			21
M	ule Index			23
In	x			25

Contents:

CONTENTS 1

2 CONTENTS

MOUSEDB CONCEPTS

Data storage for MouseDB is separated into packages which contain information about animals, and information collected about animals. There is also a separate module for timed matings of animals. This document will describe the basics of how data is stored in each of these modules.

1.1 Animal Module

Animals are tracked as individual entities, and given associations to breeding cages to follow ancestry, and strains.

1.1.1 Animal

1.1.2 Breeding Cages

A breeding cage is defined as a set of one or more male and one or more female mice. Because of this, it is not always clear who the precise parentage of an animal is. If the parentage is known, then the Mother and Father fields can be set for a particular animal.

1.1.3 Strains

1.2 Data Module

Data (or measurements) can be stored for any type of measurement. Conceptually, several pieces of data belong to an experiment (for example several mice are measured at some time) and several experiments belong to a study. Measurements can be stored independent of experiments and experiments can be performed outside of the context of a study. It is however, perfered that measurements are stored within an experiment and experiments are stored within studies as this will greatly facilitate the organization of the data.

1.2.1 Studies

In general studies are a collection of experiments. These can be grouped together on the basis of animals and/or treatment groups. A study must have at least one treatment group, which defines the animals and their conditions.

1.2.2 Experiments

An experiment is a collection of measurements for a given set of animals. In general, an experiment is defined as a number of measurements take in a given day.

1.2.3 Measurements

A measurement is an animal, an assay and a measurement value. It can be associated with an experiment, or can stand alone as an individual value. Measurements can be viewed in the context of a study, an experiment, a treatment group or an animal by going to the appropriate page.

MOUSEDB INSTALLATION

2.1 Configuration

MouseDB requires both a database and a webserver to be set up. Ideally, the database should be hosted separately from the webserver and MouseDB installation, but this is not necessary, as both can be used from the same server. If you are using a remote server for the database, it is best to set up a user for this database that can only be accessed from the webserver. If you want to set up several installations (ie for different users or different laboratories), you need separate databases and MouseDB installations for each. You will also need to set up the webserver with different addresses for each installation.

2.2 Software Dependencies

- 1. **MouseDB source code**. Download from one of the following:
- 1. http://github.com/davebridges/mousedb/downloads for the current release
- 2. http://github.com/davebridges/mousedb for the source code via Git

Downloading and/or unzipping will create a directory named mousedb. You can update to the newest revision at any time either using git or downloading and re-installing the newer version. Changing or updating software versions will not alter any saved data, but you will have to update the localsettings.py file (described below).

- 1. **Python**. Requires Version 2.6, is not yet compatible with Python 3.0. Download from http://www.python.org/download/.
- 2. Django. Download from http://www.djangoproject.com/download/
- 3. **Database software**. Typically MySQL is used, but PostgreSQL, Oracle or SQLite can also be used. You also need to install the python driver for this database (unless you are using SQLite, which is internal to Python 2.5+). See http://docs.djangoproject.com/en/dev/topics/install/database-installation Django Database Installation Documentation for more information.

2.3 Database Setup

- 1. Create a new database. You need to record the user, password, host and database name. If you are using SQLite this step is not required.
- 2. Go to localsettings_empty.py and edit the settings:
- DATABASE_ENGINE: 'mysql', 'postgresql_psycopg2' or 'sqlite3 depending on the database software used.

- DATABASE NAME: database name
- DATABASE_USER: database user
- DATABASE_PASSWORD: database password
- DATABASE_HOST: database host
- 1. Save this file as localsettings.py in the main MouseDB directory.

2.4 Web Server Setup

You need to set up a server to serve both the django installation and saved files. For the saved files. I recommend using apache for both. The preferred setup is to use Apache2 with mod_python. The following is a httpd.conf example where the code is placed in /usr/src/mousedb:

```
Alias /static /usr/src/mousedb/media
Alias /media /usr/src/mousedb/media
<Directory /usr/mousedb/media>
   Order allow, deny
   Allow from all
</Directory>
<Location "/mousedb/">
   SetHandler python-program
   PythonHandler django.core.handlers.modpython
   SetEnv DJANGO_SETTINGS_MODULE mousedb.settings
   SetEnv PYTHON_EGG_CACHE /var/www/eggs
   PythonOption django.root /mousedb
   PythonDebug On
   PythonPath "['/usr/src'] + sys.path"
   PythonInterpreter mousedb
</Location>
```

If you want to restrict access to these files, change the Allow from all directive to specific domains or ip addresses (for example Allow from 192.168.0.0/99 would allow from 192.168.0.99)

2.5 Final Configuration and User Setup

1. Go to mousedb/admin/auth/users/ and create users, selecting usernames, full names, password (or have the user set the password) and then choose group permissions.

ANIMAL DATA ENTRY

3.1 Newborn Mice or Newly Weaned Mice

- 1. Go to Breeding Cages Tab
- 2. Click on Add/Wean Pups Button
- 3. Each row is a new animal. If you accidentaly enter an extra animal, check off the delete box then submit.
- 4. Leave extra lines blank if you have less than 10 mice to enter
- 5. If you need to enter more than 10 mice, enter the first ten and submit them. Go back and enter up to 10 more animals (10 more blank spaces will appear)

3.2 Newborn Mice

- 1. Enter Breeding Cage under Cage
- 2. Enter Strain
- 3. Enter Background (normally Mixed or C57BL/6-BA unless from the LY breeding cages in which case it is C57BL/6-LY5.2)
- 4. Enter Birthdate in format YYYY-MM-DD
- 5. Enter Generation and Backcross

3.3 Weaning Mice

- 1. If not previously entered, enter data as if newborn mice
- 2. Enter gender
- 3. Enter Wean Date in format YYYY-MM-DD
- 4. Enter new Cage number for Cage

3.4 Cage Changes (Not Weaning)

1. Find mouse either from animal list or strain list

- 2. Click the edit mouse button
- 3. Change the Cage, Rack and Rack Position as Necessary

3.5 Genotyping or Ear Tagging

- 1. Find mouse either from animal list or strain list, or through breeding cage
- 2. Click the edit mouse button or the Eartag/Genotype/Cage Change/Death Button
- 3. Enter the Ear Tag and/or select the Genotype from the Pull Down List

3.6 Marking Mice as Dead

3.6.1 Dead Mice (Single Mouse)

- 1. Find mouse from animal list or strain list
- 2. Click the edit mouse button
- 3. Enter the death date in format YYYY-MM-DD
- 4. Choose Cause of Death from Pull Down List

3.6.2 Dead Mice (Several Mice)

- 1. Find mice from breeding cages
- 2. Click the Eartag/Genotype/Cage Change/Death Button
- 3. Enter the death date in format YYYY-MM-DD
- 4. Choose the Cause of Death from Pull Down List

CHAPTER

FOUR

STUDIES AND EXPERIMENTAL SETUP

Set up a new study at /mousedb/admin/data/study/ selecting animals

You must put a description and select animals in one or more treatment groups

If you have more than 2 treatment groups save the first two, then two more empty slots will appear. For animals, click on the magnifying glass then find the animal in that treatment group and click on the MouseID. The number displayed now in that field will not be the MouseID, but don't worry its just a different number to describe the mouse. To add more animals, click on the magnifying glass again and select the next animal. There should be now two numbers, separated by commas in this field. Repeat to fill all your treatment groups. You must enter a diet and environment for each treatment. The other fields are optional, and should only be used if appropriate. Ensure for pharmaceutical, you include a saline treatment group.

MEASURMENT ENTRY

5.1 Studies

If this measurement is part of a study (ie a group of experiments) then click on the plus sign beside the study field and enter in the details about the study and treatment groups. Unfortunately until i can figure out how to filter the treatment group animals in the admin interface, at each of the subsequent steps you will see all the animals in the database (soon hopefully it will only be the ones as part of the study group).

5.2 Experiment Details

- · Pick experiment date, feeding state and resarchers
- Pick animals used in this experiment (the search box will filter results)
- Fasting state, time, injections, concentration, experimentID and notes are all optional

5.3 Measurements

- There is room to enter 14 measurements. If you need more rows, enter the first 14 and select "Save and Continue Editing" and 14 more blank spots will appear.
- Each row is a measurement, so if you have glucose and weight for some animal that is two rows entered.
- For animals, click on the magnifying glass then find the animal in that treatment group and click on the MouseID. The number displayed now in that field will not be the MouseID, but don't worry its just a different number to describe the mouse.
- For values, the standard units (defined by each assay) are mg for weights, mg/dL for glucose and pg/mL for insulin). You must enter integers here (no decimal places). If you have several measurements (ie several glucose readings during a GTT, enter them all in one measurement row, separated by commas and *NO spaces*).

AUTOMATED DOCUMENTATION

6.1 Data Package

6.1.1 Models

```
class Assay (*args, **kwargs)
    Bases: django.db.models.base.Model
    Assay(id, assay, assay_slug, notes, measurement_units)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
    exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
    measurement set
class Diet (*args, **kwargs)
    Bases: django.db.models.base.Model
    Diet(id, vendor_id, description, product_id, fat_content, protein_content, carb_content, irradiated, notes)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
    exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
    treatment set
    vendor
class Environment (*args, **kwargs)
    Bases: django.db.models.base.Model
    Environment(id, building, room, temperature, humidity, notes)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
    exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
    contact
    treatment set
```

```
class Experiment (*args, **kwargs)
    Bases: django.db.models.base.Model
    Experiment(id, date, notes, experimentID, feeding_state, fasting_time, injection, concentration, study_id)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
    exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
    animals
    get_feeding_state_display(*moreargs, **morekwargs)
    get_injection_display(*moreargs, **morekwargs)
    get_next_by_date(*moreargs, **morekwargs)
    get_previous_by_date(*moreargs, **morekwargs)
    measurement_set
    researchers
    study
class Implantation (*args, **kwargs)
    Bases: django.db.models.base.Model
    Implantation(id, implant, vendor_id, product_id, notes)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
    exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
     surgeon
    treatment_set
    vendor
class Measurement (*args, **kwargs)
    Bases: django.db.models.base.Model
    Measurement(id, animal_id, experiment_id, assay_id, values)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
    exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
    animal
    assay
    experiment
class Pharmaceutical (*args, **kwargs)
    Bases: django.db.models.base.Model
    Pharmaceutical(id, drug, dose, recurrance, mode, vendor_id, notes)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
```

```
exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
    get_mode_display(*moreargs, **morekwargs)
    treatment set
    vendor
class Researcher (*args, **kwargs)
    Bases: django.db.models.base.Model
    Researcher(id, first_name, last_name, name_slug, email, active)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
    exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
    environment_set
    experiment set
    get_absolute_url (*moreargs, **morekwargs)
    implantation_set
    transplantation_set
    treatment set
class Study (*args, **kwargs)
    Bases: django.db.models.base.Model
    Study(id, description, start_date, stop_date, notes)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
    exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
    experiment_set
    strain
    treatment_set
class Transplantation (*args, **kwargs)
    Bases: django.db.models.base.Model
    Transplantation(id, tissue, transplant_date, notes)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
    exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
    donor
    get_next_by_transplant_date(*moreargs, **morekwargs)
    get_previous_by_transplant_date(*moreargs, **morekwargs)
    surgeon
```

6.1. Data Package 15

```
treatment set
class Treatment (*args, **kwargs)
    Bases: django.db.models.base.Model
    Treatment(id, treatment, study_id, diet_id, environment_id, transplantation_id, notes)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
    exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
    animals
    diet
    environment
    implantation
    pharmaceutical
    researchers
    study
    transplantation
class Vendor (*args, **kwargs)
    Bases: django.db.models.base.Model
    Vendor(id, vendor, website, email, ordering, notes)
    exception DoesNotExist
         Bases: django.core.exceptions.ObjectDoesNotExist
    exception MultipleObjectsReturned
         Bases: django.core.exceptions.MultipleObjectsReturned
    diet_set
    get_ordering_display(*moreargs, **morekwargs)
     implantation_set
    pharmaceutical_set
6.1.2 Forms
class ExperimentForm (data=None,
                                   files=None,
                                                 auto id='id %s',
                                                                   prefix=None,
                                                                                  initial=None,
                                              'django.forms.util.ErrorList'>,
                                                                                label_suffix=':',
                      error_class=<class
                      empty_permitted=False, instance=None)
    Bases: django.forms.models.ModelForm
    class Meta()
         model
             alias of Experiment
    media
```

```
files=None,
class MeasurementForm (data=None,
                                                    auto_id='id_%s',
                                                                        prefix=None,
                                                                                       initial=None,
                                                 'django.forms.util.ErrorList'>,
                        error class=<class
                                                                                     label_suffix=':',
                        empty permitted=False, instance=None)
     Bases: django.forms.models.ModelForm
     class Meta()
         model
             alias of Measurement
     media
class StudyExperimentForm (data=None, files=None,
                                                       auto_id='id_%s', prefix=None,
                                                                                       initial=None,
                             error_class=<class
                                                    'django.forms.util.ErrorList'>,
                                                                                     label_suffix=':',
                             empty_permitted=False, instance=None)
     Bases: django.forms.models.ModelForm
     class Meta()
         model
             alias of Experiment
     media
class StudyForm (data=None, files=None, auto_id='id_%s', prefix=None, initial=None, error_class=<class
                 'django.forms.util.ErrorList'>, label_suffix=':', empty_permitted=False, instance=None)
     Bases: django.forms.models.ModelForm
     class Meta()
         model
             alias of Study
     media
6.1.3 Views and URLs
add_measurement
experiment_detail
experiment_detail_all
experiment_list
study_experiment
```

6.1.4 Administrative Site Configuration

6.2 Animals Package

6.2.1 Models

This module describes the Strain, Animal, Breeding and Cage data models.

This module stores all data regarding a particular laboratory animal. Information about experimental data and timed matings are stored in the data and timed_matings packages. This module describes the database structure for each data model.

```
class Strain (*args, **kwargs)
```

Bases: django.db.models.base.Model

A data model describing a mouse strain.

This is separate from the background of a mouse. For example a ob/ob mouse on a mixed or a black-6 background still have the same strain. The background is defined in the animal and breeding cages.

```
class Animal (*args, **kwargs)
```

Bases: django.db.models.base.Model

A data model describing an animal.

This data model describes a wide variety of parameters of an experimental animal. This model is linked to the Strain and Cage models via 1:1 relationships. If the parentage of a mouse is known, this can be identified (the breeding set may not be clear on this matter). Mice are automatically marked as not alive when a Death date is provided and the object is saved.

```
class Breeding (*args, **kwargs)
```

Bases: django.db.models.base.Model

This data model stores information about a particular breeding set

A breeding set may contain one ore more males and females and must be defined via the progeny strain. For example, in the case of generating a new strain, the strain indicates the new strain not the parental strains. If the breeding set is part of a timed mating experiment, then Timed_Mating must be selected. Breeding cages are automatically inactivated upon saving when a End date is provided.

```
class Cage (*args, **kwargs)
```

Bases: django.db.models.base.Model

This data model stores information about a particular cage.

This model, which is not yet implemented will be used by both breeding and non-breeding cages and will facilitate easier tracking and storage of cages. To implement this, it will be necessary to automatically generate a new cage (if a novel barcode is entered), or to use a current cage if the barcode is already present in the database

6.2.2 Forms

```
model
            alias of Animal
    media
6.2.3 Views and URLs
animal_change
animal_detail
animal new
breeding
breeding_all
breeding_change
breeding_detail
breeding_pups
strain detail
strain_detail_all
strain_list
6.2.4 Administrative Site Configuration
class AnimalAdmin (model, admin_site)
    Bases: django.contrib.admin.options.ModelAdmin
    mark_sacrificed (request, queryset)
    media
class AnimalInline (parent_model, admin_site)
    Bases: django.contrib.admin.options.TabularInline
    model
        alias of Animal
class BreedingAdmin (model, admin_site)
    Bases: django.contrib.admin.options.ModelAdmin
    media
class CageAdmin (model, admin_site)
    Bases: django.contrib.admin.options.ModelAdmin
    media
class StrainAdmin (model, admin_site)
    Bases: django.contrib.admin.options.ModelAdmin
    media
```

class Meta()

- **6.3 Timed Mating Package**
- 6.4 Groups Package

CHAPTER

SEVEN

INDICES AND TABLES

- Index
- Module Index
- Search Page

MODULE INDEX

Α animal, 17 animal.admin, 19 animal.forms, 18 animal.models, 17 animal.urls, 19 animal.views, 19 D data, 13 data.admin, 17 ${\tt data.forms,16}$ data.models, 13 data.urls, 17 data.views, 17 G groups, 20 ${\tt timed_mating, 20}$

24 Module Index

INDEX

A	D			
add_measurement (in module data.views), 17	data (module), 13			
Animal (class in animal.models), 18	data.admin (module), 17			
animal (data.models.Measurement attribute), 14	data.forms (module), 16			
animal (module), 17	data.models (module), 13			
animal.admin (module), 19	data.urls (module), 17			
animal.forms (module), 18	data.views (module), 17			
animal.models (module), 17	Diet (class in data.models), 13			
animal.urls (module), 19	diet (data.models.Treatment attribute), 16			
animal.views (module), 19	Diet.DoesNotExist, 13			
animal_change (in module animal.views), 19	Diet.MultipleObjectsReturned, 13			
animal_detail (in module animal.views), 19	diet_set (data.models.Vendor attribute), 16			
animal_new (in module animal.views), 19	donor (data.models.Transplantation attribute), 15			
AnimalAdmin (class in animal.admin), 19	•			
AnimalChangeForm (class in animal.forms), 18	E			
AnimalChangeForm.Meta (class in animal.forms), 18	Environment (class in data.models), 13			
AnimalForm (class in animal.forms), 18	environment (data.models.Treatment attribute), 16			
AnimalForm.Meta (class in animal.forms), 18	Environment.DoesNotExist, 13			
AnimalInline (class in animal.admin), 19	Environment.MultipleObjectsReturned, 13			
animals (data.models.Experiment attribute), 14	environment_set (data.models.Researcher attribute), 15			
animals (data.models.Treatment attribute), 16	Experiment (class in data.models), 13			
Assay (class in data.models), 13	experiment (data.models.Measurement attribute), 14			
assay (data.models.Measurement attribute), 14	Experiment.DoesNotExist, 14			
Assay.DoesNotExist, 13	Experiment.MultipleObjectsReturned, 14			
Assay.MultipleObjectsReturned, 13	experiment_detail (in module data.views), 17			
5	experiment_detail_all (in module data.views), 17			
В	experiment_list (in module data.views), 17			
Breeding (class in animal.models), 18	experiment_set (data.models.Researcher attribute), 15			
breeding (in module animal.views), 19	experiment_set (data.models.Study attribute), 15			
breeding_all (in module animal.views), 19	ExperimentForm (class in data.forms), 16			
breeding_change (in module animal.views), 19	ExperimentForm.Meta (class in data.forms), 16			
breeding_detail (in module animal.views), 19				
breeding_pups (in module animal.views), 19	G			
BreedingAdmin (class in animal.admin), 19	get_absolute_url() (data.models.Researcher method), 15			
	get_feeding_state_display() (data.models.Experiment			
C	method), 14			
Cage (class in animal.models), 18	get_injection_display() (data.models.Experiment			
CageAdmin (class in animal.admin), 19	method), 14			
contact (data.models.Environment attribute), 13	get_mode_display() (data.models.Pharmaceutical method), 15			
	get next by date() (data.models.Experiment method), 14			

get_next_by_transplant_date()	R
(data.models.Transplantation method), 15	Researcher (class in data.models), 15
get_ordering_display() (data.models.Vendor method), 16	Researcher. Does Not Exist, 15
get_previous_by_date() (data.models.Experiment	Researcher.MultipleObjectsReturned, 15
method), 14	researchers (data.models.Experiment attribute), 14
get_previous_by_transplant_date()	researchers (data.models.Treatment attribute), 16
(data.models.Transplantation method), 15	6
groups (module), 20	S
I	Strain (class in animal.models), 18
Tourlandation (along in data modula) 14	strain (data.models.Study attribute), 15
Implantation (class in data.models), 14 implantation (data.models.Treatment attribute), 16	strain_detail (in module animal.views), 19
Implantation.DoesNotExist, 14	strain_detail_all (in module animal.views), 19
Implantation.MultipleObjectsReturned, 14	strain_list (in module animal.views), 19
implantation_set (data.models.Researcher attribute), 15	StrainAdmin (class in animal.admin), 19 Study (class in data.models), 15
implantation_set (data.models.Vendor attribute), 16	study (data.models, Experiment attribute), 14
	study (data.models.Experiment attribute), 14
M	Study. Does Not Exist, 15
mark_sacrificed() (animal.admin.AnimalAdmin method),	Study.MultipleObjectsReturned, 15
19	study_experiment (in module data.views), 17
Measurement (class in data.models), 14	StudyExperimentForm (class in data.forms), 17
Measurement.DoesNotExist, 14	StudyExperimentForm.Meta (class in data.forms), 17
Measurement.MultipleObjectsReturned, 14	StudyForm (class in data.forms), 17
measurement_set (data.models.Assay attribute), 13	StudyForm.Meta (class in data.forms), 17
measurement_set (data.models.Experiment attribute), 14	surgeon (data.models.Implantation attribute), 14
MeasurementForm (class in data forms), 16	surgeon (data.models.Transplantation attribute), 15
MeasurementForm.Meta (class in data.forms), 17 media (animal.admin.AnimalAdmin attribute), 19	T
media (animal.admin.BreedingAdmin attribute), 19	
media (animal.admin.CageAdmin attribute), 19	timed_mating (module), 20 Transplantation (along in data module), 15
media (animal.admin.StrainAdmin attribute), 19	Transplantation (class in data.models), 15 transplantation (data.models.Treatment attribute), 16
media (animal.forms.AnimalChangeForm attribute), 18	Transplantation. Does Not Exist, 15
media (animal.forms.AnimalForm attribute), 19	Transplantation.MultipleObjectsReturned, 15
media (data.forms.ExperimentForm attribute), 16	transplantation_set (data.models.Researcher attribute), 15
media (data.forms.MeasurementForm attribute), 17	Treatment (class in data.models), 16
media (data.forms.StudyExperimentForm attribute), 17	Treatment.DoesNotExist, 16
media (data.forms.StudyForm attribute), 17	Treatment.MultipleObjectsReturned, 16
model (animal.admin.AnimalInline attribute), 19	treatment_set (data.models.Diet attribute), 13
model (animal.forms.AnimalChangeForm.Meta attribute), 18	treatment_set (data.models.Environment attribute), 13
model (animal.forms.AnimalForm.Meta attribute), 19	treatment_set (data.models.Implantation attribute), 14
model (data.forms.ExperimentForm.Meta attribute), 16	treatment_set (data.models.Pharmaceutical attribute), 15
model (data.forms.MeasurementForm.Meta attribute), 17	treatment_set (data.models.Researcher attribute), 15 treatment_set (data.models.Study attribute), 15
model (data.forms.StudyExperimentForm.Meta at-	treatment_set (data.models.Transplantation attribute), 15
tribute), 17	treatment_set (data.models.fransplantation attribute), 13
model (data.forms.StudyForm.Meta attribute), 17	V
П	Vendor (class in data.models), 16
P	vendor (data.models.Diet attribute), 13
Pharmaceutical (class in data.models), 14	vendor (data.models.Implantation attribute), 14
pharmaceutical (data.models.Treatment attribute), 16	vendor (data.models.Pharmaceutical attribute), 15
Pharmaceutical.DoesNotExist, 14	Vendor.DoesNotExist, 16
Pharmaceutical.MultipleObjectsReturned, 14	Vendor.MultipleObjectsReturned, 16
pharmaceutical_set (data.models.Vendor attribute), 16	

26 Index