
MouseDB Documentation

Release 0.1

Dave Bridges, Ph.D.

May 16, 2010

CONTENTS

1	MouseDB Concepts	3
1.1	Animal Module	3
1.2	Data Module	4
1.3	Timed Matings Module	4
1.4	Groups Module	4
2	MouseDB Installation	5
2.1	Configuration	5
2.2	Software Dependencies	5
2.3	Database Setup	5
2.4	Web Server Setup	6
2.5	Final Configuration and User Setup	6
3	Animal Data Entry	7
3.1	Newborn Mice or Newly Weaned Mice	7
3.2	Newborn Mice	7
3.3	Weaning Mice	7
3.4	Cage Changes (Not Weaning)	7
3.5	Genotyping or Ear Tagging	8
3.6	Marking Mice as Dead	8
4	Studies and Experimental Setup	9
5	Measurment Entry	11
5.1	Studies	11
5.2	Experiment Details	11
5.3	Measurements	11
6	Automated Documentation	13
6.1	Data Package	13
6.2	Animals Package	18
6.3	Timed Mating Package	22
6.4	Groups Package	24
7	Indices and tables	27
	Module Index	29
	Index	31

Contents:

MOUSEDB CONCEPTS

Data storage for MouseDB is separated into packages which contain information about animals, and information collected about animals. There is also a separate module for timed matings of animals. This document will describe the basics of how data is stored in each of these modules.

1.1 Animal Module

Animals are tracked as individual entities, and given associations to breeding cages to follow ancestry, and strains.

1.1.1 Animal

Most parameters about an animal are set within the animal object. Here is where the animals strain, breeding, parentage and many other parameters are included. Animals have foreignkey relationships with both Strain and Breeding, so an animal may only belong to one of each of those. As an example, a mouse cannot come from more than one Breeding set, and cannot belong to more than one strain.

Backcrosses and Generations

For this software, optional tracking of backcrosses and generations is available and is stored as an attribute of an animal. When an inbred cross is made against a pure background, the backcross increases by 1. When a heterozygote cross is made, the generation increases by one. As an example, for every time a mouse in a C57/BL6 background is crossed against a wildtype C57/B6 mouse, the backcross (but not the generation) increases by one. For every time a mutant strain is crosses against itself (either vs a heterozygote or homozygote of that strain), the generation will increase by one. Backcrosses should typically be performed against a separate colony of purebred mouse, rather than against wild-type alleles of the mutant strain.

1.1.2 Breeding Cages

A breeding cage is defined as a set of one or more male and one or more female mice. Because of this, it is not always clear who the precise parentage of an animal is. If the parentage is known, then the Mother and Father fields can be set for a particular animal. In the case of Active, if an End field is specified, then the Active field is set to False. In the case of Cage, if a Cage is provided, and animals are specified under Male or Females for a Breeding object, then the Cage field for those animals is set to that of the breeding cage. The same is true for both Rack and Rack Position.

1.1.3 Strains

A strain is a set of mice with a similar genetics. Importantly strains are separated from Backgrounds. For example, one might have mice with the genotype ob/ob but these mice may be in either a C57-Black6 or a mixed background. This difference is set at the individual animal level. The result of this is that a query for a particular strain may then need to be filtered to a specific background.

1.2 Data Module

Data (or measurements) can be stored for any type of measurement. Conceptually, several pieces of data belong to an experiment (for example several mice are measured at some time) and several experiments belong to a study. Measurements can be stored independent of experiments and experiments can be performed outside of the context of a study. It is however, preferred that measurements are stored within an experiment and experiments are stored within studies as this will greatly facilitate the organization of the data.

1.2.1 Studies

In general studies are a collection of experiments. These can be grouped together on the basis of animals and/or treatment groups. A study must have at least one treatment group, which defines the animals and their conditions.

1.2.2 Experiments

An experiment is a collection of measurements for a given set of animals. In general, an experiment is defined as a number of measurements taken in a given day.

1.2.3 Measurements

A measurement is an animal, an assay and a measurement value. It can be associated with an experiment, or can stand alone as an individual value. Measurements can be viewed in the context of a study, an experiment, a treatment group or an animal by going to the appropriate page.

1.3 Timed Matings Module

Timed matings are a specific type of breeding set. Generally, for these experiments a mating cage is set up and pregnancy is defined by a plug event. Based on this information, the age of an embryo can be estimated. When a breeding cage is defined, one option is to set this cage as a timed mating cage (ie `Timed_Mating=True`). If this is the case, then a plug event can be registered and recorded for this mating set. If the mother gives birth then this cage is implicitly set as a normal breeding cage.

1.4 Groups Module

This app defines generic Group and License information for a particular installation of MouseDB. Because every page on this site identifies both the Group and data restrictions, at a minimum, group information must be provided upon installation (see installation instructions).

MOUSEDDB INSTALLATION

2.1 Configuration

MouseDB requires both a database and a webserver to be set up. Ideally, the database should be hosted separately from the webserver and MouseDB installation, but this is not necessary, as both can be used from the same server. If you are using a remote server for the database, it is best to set up a user for this database that can only be accessed from the webserver. If you want to set up several installations (ie for different users or different laboratories), you need separate databases and MouseDB installations for each. You will also need to set up the webserver with different addresses for each installation.

2.2 Software Dependencies

1. **MouseDB source code.** Download from one of the following:
 1. <http://github.com/davebridges/mousedb/downloads> for the current release
 2. <http://github.com/davebridges/mousedb> for the source code via Git

Downloading and/or unzipping will create a directory named mousedb. You can update to the newest revision at any time either using git or downloading and re-installing the newer version. Changing or updating software versions will not alter any saved data, but you will have to update the localsettings.py file (described below).

1. **Python.** Requires Version 2.6, is not yet compatible with Python 3.0. Download from <http://www.python.org/download/>.
2. **Django.** Download from <http://www.djangoproject.com/download/>
3. **Database software.** Typically MySQL is used, but PostgreSQL, Oracle or SQLite can also be used. You also need to install the python driver for this database (unless you are using SQLite, which is internal to Python 2.5+). See <http://docs.djangoproject.com/en/dev/topics/install/database-installation> - Django Database Installation Documentation for more information.

2.3 Database Setup


1. Create a new database. You need to record the user, password, host and database name. If you are using SQLite this step is not required.
2. Go to localsettings_empty.py and edit the settings:
 - DATABASE_ENGINE: 'mysql', 'postgresql_psycopg2' or 'sqlite3' depending on the database software used.

- DATABASE_NAME: database name
- DATABASE_USER: database user
- DATABASE_PASSWORD: database password
- DATABASE_HOST: database host

1. Save this file as localsettings.py in the main MouseDB directory.

2.4 Web Server Setup

You need to set up a server to serve both the django installation and saved files. For the saved files, I recommend using apache for both. The preferred setup is to use Apache2 with mod_python. The following is a httpd.conf example where the code is placed in /usr/src/mousedb:



If you want to restrict access to these files, change the Allow from all directive to specific domains or ip addresses (for example Allow from 192.168.0.0/99 would allow from 192.168.0.0 to 192.168.0.99)

2.5 Final Configuration and User Setup

Go to a command prompt, navigate to inside the mousedb directory and enter the following to get to a python prompt:



From the command prompt enter the following where **groupname** is the name of your research group:

```
from mousedb.groups.models import  
    =  
    ="groupname"
```

Go to mousedb/admin/auth/users/ and create users, selecting usernames, full names, password (or have the user set the password) and then choose group permissions.

ANIMAL DATA ENTRY

3.1 Newborn Mice or Newly Weaned Mice

1. Go to Breeding Cages Tab
2. Click on Add/Wean Pups Button
3. Each row is a new animal. If you accidentally enter an extra animal, check off the delete box then submit.
4. Leave extra lines blank if you have less than 10 mice to enter
5. If you need to enter more than 10 mice, enter the first ten and submit them. Go back and enter up to 10 more animals (10 more blank spaces will appear)

3.2 Newborn Mice

1. Enter Breeding Cage under Cage
2. Enter Strain
3. Enter Background (normally Mixed or C57BL/6-BA unless from the LY breeding cages in which case it is C57BL/6-LY5.2)
4. Enter Birthdate in format YYYY-MM-DD
5. Enter Generation and Backcross

3.3 Weaning Mice

1. If not previously entered, enter data as if newborn mice
2. Enter gender
3. Enter Wean Date in format YYYY-MM-DD
4. Enter new Cage number for Cage

3.4 Cage Changes (Not Weaning)

1. Find mouse either from animal list or strain list

2. Click the edit mouse button
3. Change the Cage, Rack and Rack Position as Necessary

3.5 Genotyping or Ear Tagging

1. Find mouse either from animal list or strain list, or through breeding cage
2. Click the edit mouse button or the Eartag/Genotype/Cage Change/Death Button
3. Enter the Ear Tag and/or select the Genotype from the Pull Down List

3.6 Marking Mice as Dead

3.6.1 Dead Mice (Single Mouse)

1. Find mouse from animal list or strain list
2. Click the edit mouse button
3. Enter the death date in format YYYY-MM-DD
4. Choose Cause of Death from Pull Down List

3.6.2 Dead Mice (Several Mice)

1. Find mice from breeding cages
2. Click the Eartag/Genotype/Cage Change/Death Button
3. Enter the death date in format YYYY-MM-DD
4. Choose the Cause of Death from Pull Down List

STUDIES AND EXPERIMENTAL SETUP

Set up a new study at </mousedb/admin/data/study/> selecting animals

You must put a description and select animals in one or more treatment groups

If you have more than 2 treatment groups save the first two, then two more empty slots will appear. For animals, click on the magnifying glass then find the animal in that treatment group and click on the MouseID. The number displayed now in that field will not be the MouseID, but don't worry its just a different number to describe the mouse. To add more animals, click on the magnifying glass again and select the next animal. There should be now two numbers, separated by commas in this field. Repeat to fill all your treatment groups. You must enter a diet and environment for each treatment. The other fields are optional, and should only be used if appropriate. Ensure for pharmaceutical, you include a saline treatment group.

MEASUREMENT ENTRY

5.1 Studies

If this measurement is part of a study (ie a group of experiments) then click on the plus sign beside the study field and enter in the details about the study and treatment groups. Unfortunately until i can figure out how to filter the treatment group animals in the admin interface, at each of the subsequent steps you will see all the animals in the database (soon hopefully it will only be the ones as part of the study group).

5.2 Experiment Details

- Pick experiment date, feeding state and resarchers
- Pick animals used in this experiment (the search box will filter results)
- Fasting state, time, injections, concentration, experimentID and notes are all optional

5.3 Measurements

- There is room to enter 14 measurements. If you need more rows, enter the first 14 and select “Save and Continue Editing” and 14 more blank spots will appear.
- Each row is a measurement, so if you have glucose and weight for some animal that is two rows entered.
- For animals, click on the magnifying glass then find the animal in that treatment group and click on the MouseID. The number displayed now in that field will not be the MouseID, but don’t worry its just a different number to describe the mouse.
- For values, the standard units (defined by each assay) are mg for weights, mg/dL for glucose and pg/mL for insulin). You must enter integers here (no decimal places). If you have several measurements (ie several glucose readings during a GTT, enter them all in one measurement row, separated by commas and *NO spaces*).

AUTOMATED DOCUMENTATION

6.1 Data Package

The data module describes the conditions and collection of data regarding experimental animals.

Data (or measurements) can be stored for any type of measurement. Conceptually, several pieces of data belong to an experiment (for example several mice are measured at some time) and several experiments belong to a study. Measurements can be stored independent of experiments and experiments can be performed outside of the context of a study. It is however, preferred that measurements are stored within an experiment and experiments are stored within studies as this will greatly facilitate the organization of the data.

6.1.1 Studies

In general studies are a collection of experiments. These can be grouped together on the basis of animals and/or treatment groups. A study must have at least one treatment group, which defines the animals and their conditions.

6.1.2 Experiments

An experiment is a collection of measurements for a given set of animals. In general, an experiment is defined as a number of measurements taken in a given day.

6.1.3 Measurements

A measurement is an animal, an assay and a measurement value. It can be associated with an experiment, or can stand alone as an individual value. Measurements can be viewed in the context of a study, an experiment, a treatment group or an animal by going to the appropriate page.

6.1.4 Models

```
class Assay (*args, **kwargs)
    Bases: django.db.models.base.Model

    Assay(id, assay, assay_slug, notes, measurement_units)

    exception DoesNotExist
        Bases: django.core.exceptions.ObjectDoesNotExist

    exception MultipleObjectsReturned
        Bases: django.core.exceptions.MultipleObjectsReturned
```

```

    measurement_set

class Diet (*args, **kwargs)
    Bases: django.db.models.base.Model

    Diet(id, vendor_id, description, product_id, fat_content, protein_content, carb_content, irradiated, notes)

    exception DoesNotExist
        Bases: django.core.exceptions.ObjectDoesNotExist

    exception MultipleObjectsReturned
        Bases: django.core.exceptions.MultipleObjectsReturned

    treatment_set

    vendor

class Environment (*args, **kwargs)
    Bases: django.db.models.base.Model

    Environment(id, building, room, temperature, humidity, notes)

    exception DoesNotExist
        Bases: django.core.exceptions.ObjectDoesNotExist

    exception MultipleObjectsReturned
        Bases: django.core.exceptions.MultipleObjectsReturned

    contact

    treatment_set

class Experiment (*args, **kwargs)
    Bases: django.db.models.base.Model

    Experiment(id, date, notes, experimentID, feeding_state, fasting_time, injection, concentration, study_id)

    exception DoesNotExist
        Bases: django.core.exceptions.ObjectDoesNotExist

    exception MultipleObjectsReturned
        Bases: django.core.exceptions.MultipleObjectsReturned

    animals

    get_absolute_url (*moreargs, **morekwargs)

    get_feeding_state_display (*moreargs, **morekwargs)

    get_injection_display (*moreargs, **morekwargs)

    get_next_by_date (*moreargs, **morekwargs)

    get_previous_by_date (*moreargs, **morekwargs)

    measurement_set

    researchers

    study

class Implantation (*args, **kwargs)
    Bases: django.db.models.base.Model

    Implantation(id, implant, vendor_id, product_id, notes)

    exception DoesNotExist
        Bases: django.core.exceptions.ObjectDoesNotExist

```

```

exception MultipleObjectsReturned
    Bases: django.core.exceptions.MultipleObjectsReturned

surgeon

treatment_set

vendor

class Measurement (*args, **kwargs)
    Bases: django.db.models.base.Model
    Measurement(id, animal_id, experiment_id, assay_id, values)

exception DoesNotExist
    Bases: django.core.exceptions.ObjectDoesNotExist

exception MultipleObjectsReturned
    Bases: django.core.exceptions.MultipleObjectsReturned

animal

assay

experiment

class Pharmaceutical (*args, **kwargs)
    Bases: django.db.models.base.Model
    Pharmaceutical(id, drug, dose, recurrence, mode, vendor_id, notes)

exception DoesNotExist
    Bases: django.core.exceptions.ObjectDoesNotExist

exception MultipleObjectsReturned
    Bases: django.core.exceptions.MultipleObjectsReturned

get_mode_display (*moreargs, **morekwargs)

treatment_set

vendor

class Researcher (*args, **kwargs)
    Bases: django.db.models.base.Model
    Researcher(id, first_name, last_name, name_slug, email, active)

exception DoesNotExist
    Bases: django.core.exceptions.ObjectDoesNotExist

exception MultipleObjectsReturned
    Bases: django.core.exceptions.MultipleObjectsReturned

environment_set

experiment_set

get_absolute_url (*moreargs, **morekwargs)

implantation_set

transplantation_set

treatment_set

```

```

class Study (*args, **kwargs)
    Bases: django.db.models.base.Model

    Study(id, description, start_date, stop_date, notes)

    exception DoesNotExist
        Bases: django.core.exceptions.ObjectDoesNotExist

    exception MultipleObjectsReturned
        Bases: django.core.exceptions.MultipleObjectsReturned

    experiment_set

    get_absolute_url (*moreargs, **morekwargs)

    strain

    treatment_set

class Transplantation (*args, **kwargs)
    Bases: django.db.models.base.Model

    Transplantation(id, tissue, transplant_date, notes)

    exception DoesNotExist
        Bases: django.core.exceptions.ObjectDoesNotExist

    exception MultipleObjectsReturned
        Bases: django.core.exceptions.MultipleObjectsReturned

    donor

    get_next_by_transplant_date (*moreargs, **morekwargs)

    get_previous_by_transplant_date (*moreargs, **morekwargs)

    surgeon

    treatment_set

class Treatment (*args, **kwargs)
    Bases: django.db.models.base.Model

    Treatment(id, treatment, study_id, diet_id, environment_id, transplantation_id, notes)

    exception DoesNotExist
        Bases: django.core.exceptions.ObjectDoesNotExist

    exception MultipleObjectsReturned
        Bases: django.core.exceptions.MultipleObjectsReturned

    animals

    diet

    environment

    get_absolute_url (*moreargs, **morekwargs)

    implantation

    pharmaceutical

    researchers

    study

    transplantation

```

```

class Vendor (*args, **kwargs)
    Bases: django.db.models.base.Model

    Vendor(id, vendor, website, email, ordering, notes)

    exception DoesNotExist
        Bases: django.core.exceptions.ObjectDoesNotExist

    exception MultipleObjectsReturned
        Bases: django.core.exceptions.MultipleObjectsReturned

    diet_set

    get_ordering_display (*moreargs, **morekwargs)

    implantation_set

    pharmaceutical_set

```

6.1.5 Forms

```

class ExperimentForm (data=None, files=None, auto_id='id_%s', prefix=None, initial=None,
                      error_class=<class 'django.forms.util.ErrorList'>, label_suffix=':',
                      empty_permitted=False, instance=None)
    Bases: django.forms.models.ModelForm

```

```

class Meta ()

```

```

    model
        alias of Experiment

```

```

media

```

```

class MeasurementForm (data=None, files=None, auto_id='id_%s', prefix=None, initial=None,
                       error_class=<class 'django.forms.util.ErrorList'>, label_suffix=':',
                       empty_permitted=False, instance=None)
    Bases: django.forms.models.ModelForm

```

Form definition for adding and editing measurements.

This form excludes experiment, which must be passed as a filtering parameter from the view. This form is used for formsets to add or modify measurements from within an experiment.

```

class Meta ()

```

```

    model
        alias of Measurement

```

```

media

```

```

class StudyExperimentForm (data=None, files=None, auto_id='id_%s', prefix=None, initial=None,
                           error_class=<class 'django.forms.util.ErrorList'>, label_suffix=':',
                           empty_permitted=False, instance=None)
    Bases: django.forms.models.ModelForm

```

```

class Meta ()

```

```

    model
        alias of Experiment

```

```

media

```

```
class StudyForm (data=None, files=None, auto_id='id_%s', prefix=None, initial=None, error_class=<class  
    'django.forms.util.ErrorList'>, label_suffix=':', empty_permitted=False, instance=None)  
    Bases: django.forms.models.ModelForm
```

```
    class Meta ()
```

```
        model
```

```
            alias of Study
```

```
        media
```

```
class TreatmentForm (data=None, files=None, auto_id='id_%s', prefix=None, initial=None, error_class=<class  
    'django.forms.util.ErrorList'>, label_suffix=':', empty_permitted=False, instance=None)  
    Bases: django.forms.models.ModelForm
```

Form class for study treatment groups.

In the case of studies, animals are defined in the treatment group rather than in the study group. A treatment consists of a study, a set of animals and the conditions which define that treatment. This includes related fields for environment, diet, implants and transplants.

```
    class Meta ()
```

```
        model
```

```
            alias of Treatment
```

```
        media
```

6.1.6 Views and URLs

```
add_measurement
```

```
experiment_detail
```

```
experiment_detail_all
```

```
experiment_list
```

```
study_experiment
```

6.1.7 Administrative Site Configuration

6.2 Animals Package

The animal app contains and controls the display of data about animals.

Animals are tracked as individual entities, and given associations to breeding cages to follow ancestry, and strains.

6.2.1 Animal

Most parameters about an animal are set within the animal object. Here is where the animals strain, breeding, parentage and many other parameters are included. Animals have foreignkey relationships with both Strain and Breeding, so an animal may only belong to one of each of those. As an example, a mouse cannot come from more than one Breeding set, and cannot belong to more than one strain.

Backcrosses and Generations

For this software, optional tracking of backcrosses and generations is available and is stored as an attribute of an animal. When an inbred cross is made against a pure background, the backcross increases by 1. When a heterozygote cross is made, the generation increases by one. As an example, for every time a mouse in a C57/BL6 background is crossed against a wildtype C57/B6 mouse, the backcross (but not the generation) increases by one. For every time a mutant strain is crosses against itself (either vs a heterozygote or homozygote of that strain), the generation will increase by one. Backcrosses should typically be performed against a separate colony of purebred mouse, rather than against wild-type alleles of the mutant strain.

6.2.2 Breeding Cages

A breeding cage is defined as a set of one or more male and one or more female mice. Because of this, it is not always clear who the precise parentage of an animal is. If the parentage is known, then the Mother and Father fields can be set for a particular animal.

6.2.3 Strains

A strain is a set of mice with a similar genetics. Importantly strains are separated from Backgrounds. For example, one might have mice with the genotype ob/ob but these mice may be in either a C57-Black6 or a mixed background. This difference is set at the individual animal level. The result of this is that a query for a particular strain may then need to be filtered to a specific background.

6.2.4 Models

This module describes the Strain, Animal, Breeding and Cage data models.

This module stores all data regarding a particular laboratory animal. Information about experimental data and timed matings are stored in the data and timed_matings packages. This module describes the database structure for each data model.

```
class Strain (*args, **kwargs)
```

```
    Bases: django.db.models.base.Model
```

A data model describing a mouse strain.

This is separate from the background of a mouse. For example a ob/ob mouse on a mixed or a black-6 background still have the same strain. The background is defined in the animal and breeding cages. Strain and Strain_slug are required.

```
class Animal (*args, **kwargs)
```

```
    Bases: django.db.models.base.Model
```

A data model describing an animal.

This data model describes a wide variety of parameters of an experimental animal. This model is linked to the Strain and Cage models via 1:1 relationships. If the parentage of a mouse is known, this can be identified (the breeding set may not be clear on this matter). Mice are automatically marked as not alive when a Death date is provided and the object is saved. Strain, Background and Genotype are required field. By default, querysets are ordered first by strain then by MouseID.

```
class Breeding (*args, **kwargs)
```

```
    Bases: django.db.models.base.Model
```

This data model stores information about a particular breeding set

A breeding set may contain one or more males and females and must be defined via the progeny strain. For example, in the case of generating a new strain, the strain indicates the new strain not the parental strains. A breeding cage is defined as one male with one or more females. If the breeding set is part of a timed mating experiment, then `Timed_Mating` must be selected. Breeding cages are automatically inactivated upon saving when an End date is provided. The only required field is `Strain`. By default, querysets are ordered by `Strain`, then `Start`.

```
class Cage (*args, **kwargs)
    Bases: django.db.models.base.Model
```

This data model stores information about a particular cage.

This model, which is not yet implemented will be used by both breeding and non-breeding cages and will facilitate easier tracking and storage of cages. To implement this, it will be necessary to automatically generate a new cage (if a novel barcode is entered), or to use a current cage if the barcode is already present in the database

6.2.5 Forms

Forms for use in manipulating objects in the animal app.

```
class AnimalForm (data=None, files=None, auto_id='id_%s', prefix=None, initial=None, error_class=<class
    'django.forms.util.ErrorList'>, label_suffix=':', empty_permitted=False, instance=None)
    Bases: django.forms.models.ModelForm
```

This modelform provides fields for modifying animal data.

It includes all fields except `CageID` (will be deprecated) and `Alive` (which is automatically set upon saving the animal). This form also automatically loads javascript and css for the datepicker jquery-ui widget.

```
class Media ()
```

```
class Meta ()
```

```
    model
        alias of Animal
```

```
    media
```

```
class BreedingForm (data=None, files=None, auto_id='id_%s', prefix=None, initial=None, error_class=<class
    'django.forms.util.ErrorList'>, label_suffix=':', empty_permitted=False, instance=None)
    Bases: django.forms.models.ModelForm
```

This form provides most fields for creating and entering breeding cage data.

This form is used from the url `/mousedb/breeding/new` and is a generic create view. The only excluded field is `CageID`, as this feature will be deprecated.

```
class Meta ()
```

```
    model
        alias of Breeding
```

```
    media
```

6.2.6 Views and URLs

These views define template redirects for the animal app.

This module contains only non-generic views. Several generic views are also used and are defined in `animal/urls/`.

animal_detail

This view displays specific details about an animal.

It takes a request in the form `animal/(id)/`, `mice/(id)` or `mouse/(id)/` and renders the detail page for that mouse. The request is defined for `id` not `MouseID` (or `barcode`) because this allows for details to be displayed for mice without barcode identification. Therefore care must be taken that `animal/4932` is `id=4932` and not `barcode=4932`. The animal name is defined at the top of the page. This page is restricted to logged-in users.

breeding

This view displays a list of breeding sets.

It takes a request in the form `/breeding/` and lists all currently active breeding sets. This page is restricted to logged-in users.

breeding_all

This view displays a list of breeding sets.

It takes a request in the form `/breeding/all` and lists all breeding sets (both active and inactive). This page is restricted to logged-in users.

breeding_change

This view is used to generate a form by which to change pups which belong to a particular breeding set.

This view typically is used to modify existing pups. This might include marking animals as sacrificed, entering genotype or marking information or entering movement of mice to another cage. It is used to show and modify several animals at once. It takes a request in the form `/breeding/(breeding_id)/change/` and returns a form specific to the breeding set defined in `breeding_id`. `breeding_id` is the background identification number of the breeding set and does not refer to the barcode of any breeding cage. This view returns a formset in which one row represents one animal. To add extra animals to a breeding set use `/breeding/(breeding_id)/pups/`. This view is restricted to those with the permission `animal.change_animal`.

breeding_detail

This view displays specific details about a breeding set.

It takes a request in the form `/breeding/(breeding_id)/all` and renders the detail page for that breeding set. The `breeding_id` refers to the background id of the breeding set, and not the breeding cage barcode. This view also passes along a dictionary of all pups belonging to that breeding set. This page is restricted to logged-in users.

breeding_pups

This view is used to generate a form by which to add pups which belong to a particular breeding set.

This view is intended to be used to add initial information about pups, including eartag, genotype, gender and birth or weaning information. It takes a request in the form `/breeding/(breeding_id)/pups/` and returns a form specific to the breeding set defined in `breeding_id`. `breeding_id` is the background identification number of the breeding set and does not refer to the barcode of any breeding cage. This view is restricted to those with the permission `animal.add_animal`.

strain_detail

This view displays specific details about a strain.

It takes a request in the form `strain/(strain_slug)/` and renders the detail page for that strain. This view also passes along a dictionary of alive animals belonging to that strain. This page is restricted to logged-in users.

strain_detail_all

This view displays specific details about a strain.

It takes a request in the form `/strain/(strain_slug)/all` and renders the detail page for that strain. This view also passes along a dictionary of all animals belonging to that strain. This page is restricted to logged-in users.

strain_list

This view presents a list of strains currently present in the database and annotates this list with a count of alive and total animals.

This view redirects from a /strain/ request. This view is restricted to logged-in users.

6.2.7 Administrative Site Configuration

Admin site settings for the animal app.

class AnimalAdmin (*model, admin_site*)

Bases: `django.contrib.admin.options.ModelAdmin`

Provides parameters for animal objects within the admin interface.

mark_sacrificed (*request, queryset*)

An admin action for marking several animals as sacrificed.

This action sets the selected animals as Alive=False, Death=today and Cause_of_Death as sacrificed. To use other paramters, mice muse be individually marked as sacrificed. This admin action also shows as the output the number of mice sacrificed.

media

class AnimalInline (*parent_model, admin_site*)

Bases: `django.contrib.admin.options.TabularInline`

Provides an inline tabular formset for animal objects.

Currently used with the breeding admin page.

model

alias of `Animal`

class BreedingAdmin (*model, admin_site*)

Bases: `django.contrib.admin.options.ModelAdmin`

Settings in the admin interface for dealing with Breeding objects.

This interface also includes an form for adding objects associated with this breeding cage.

mark_deactivated (*request, queryset*)

An admin action for marking several cages as inactive.

This action sets the selected cages as Active=False and Death=today. This admin action also shows as the output the number of mice sacrificed.

media

class StrainAdmin (*model, admin_site*)

Bases: `django.contrib.admin.options.ModelAdmin`

Settings in the admin interface for dealing with Strain objects.

media

6.3 Timed Mating Package

This package defines the `timed_mating` app.

Timed matings are a specific type of breeding set. Generally, for these experiments a mating cage is set up and pregnancy is defined by a plug event. Based on this information, the age of an embryo can be estimated. When a breeding cage is defined, one option is to set this cage as a timed mating cage (ie `Timed_Mating=True`). If this is the case, then a plug event can be registered and recorded for this mating set. If the mother gives birth then this cage is implicitly set as a normal breeding cage.

6.3.1 Models

This defines the data model for the `timed_mating` app.

Currently the only data model is for `PlugEvents`.

class `PlugEvents` (**args*, ***kwargs*)

Bases: `django.db.models.base.Model`

This defines the model for `PlugEvents`.

A `PlugEvent` requires a date. All other fields are optional. Upon observation of a plug event, the `PlugDate`, `Breeding Cage`, `Femalem`, `Male`, `Researcher` and `Notes` can be set. Upon sacrifice of the mother, then genotyped alive and dead embryos can be entered, along with the `SacrificeDate`, `Researcher` and `Notes`.

Breeding

exception `DoesNotExist`

Bases: `django.core.exceptions.ObjectDoesNotExist`

exception `MultipleObjectsReturned`

Bases: `django.core.exceptions.MultipleObjectsReturned`

`PlugFemale`

`PlugMale`

`Researcher`

`get_absolute_url` (**moreargs*, ***morekwargs*)

`get_next_by_PlugDate` (**moreargs*, ***morekwargs*)

`get_previous_by_PlugDate` (**moreargs*, ***morekwargs*)

`save` ()

Over-rides the default save function for `PlugEvents`.

If a sacrifice date is set for an object in this model, then `Active` is set to `False`.

class `PlugEvents` (**args*, ***kwargs*)

Bases: `django.db.models.base.Model`

This defines the model for `PlugEvents`.

A `PlugEvent` requires a date. All other fields are optional. Upon observation of a plug event, the `PlugDate`, `Breeding Cage`, `Femalem`, `Male`, `Researcher` and `Notes` can be set. Upon sacrifice of the mother, then genotyped alive and dead embryos can be entered, along with the `SacrificeDate`, `Researcher` and `Notes`.

6.3.2 Forms

This package describes forms used by the `Timed Mating` app.

```
class BreedingPlugForm (data=None, files=None, auto_id='id_%s', prefix=None, initial=None,
                        error_class=<class 'django.forms.util.ErrorList'>, label_suffix=':',
                        empty_permitted=False, instance=None)
Bases: django.forms.models.ModelForm
```

This form is used to enter Plug Events from a specific breeding cage.

```
class Meta ()
```

```
    model
        alias of PlugEvents
```

```
media
```

6.3.3 Views and URLs

This package defines custom views for the timed_mating application.

Currently all views are generic CRUD views except for the view in which a plug event is defined from a breeding cage.

breeding_plugevent

This view defines a form for adding new plug events from a breeding cage.

This form requires a breeding_id from a breeding set and restricts the PlugFemale and PlugMale to animals that are defined in that breeding cage.

change_plugevents

create_plugevents

delete_plugevents

limited_object_detail

limited_object_list

6.3.4 Administrative Site Configuration

Settings to control the admin interface for the timed_mating app.

This file defines a PlugEventsAdmin object to enter parameters about individual plug events/

```
class PlugEventsAdmin (model, admin_site)
```

```
    Bases: django.contrib.admin.options.ModelAdmin
```

This class defines the admin interface for the PlugEvents model.

```
media
```

6.4 Groups Package

This package defines the Group application. This app defines generic Group and License information for a particular installation of MouseDB. Because every page on this site identifies both the Group and data restrictions, at a minimum, group information must be provided upon installation (see installation instructions).

6.4.1 Models

class Group (*args, **kwargs)

Bases: django.db.models.base.Model

This defines the data structure for the Group model.

The only required field is group. All other fields (group_slug, group_url, license, contact_title, contact_first, contact_last and contact_email) are optional.

exception DoesNotExist

Bases: django.core.exceptions.ObjectDoesNotExist

exception MultipleObjectsReturned

Bases: django.core.exceptions.MultipleObjectsReturned

get_contact_title_display (*moreargs, **morekwargs)

license

class License (*args, **kwargs)

Bases: django.db.models.base.Model

This defines the data structure for the License model.

The only required field is license. If the contents of this installation are being made available using some licencing criteria this can either be defined in the notes field, or in an external website.

exception DoesNotExist

Bases: django.core.exceptions.ObjectDoesNotExist

exception MultipleObjectsReturned

Bases: django.core.exceptions.MultipleObjectsReturned

group_set

class Group (*args, **kwargs)

Bases: django.db.models.base.Model

This defines the data structure for the Group model.

The only required field is group. All other fields (group_slug, group_url, license, contact_title, contact_first, contact_last and contact_email) are optional.

6.4.2 Views and URLs

6.4.3 Administrative Site Configuration

class GroupAdmin (model, admin_site)

Bases: django.contrib.admin.options.ModelAdmin

Defines the admin interface for Groups.

Currently set as default.

media

class LicenseAdmin (model, admin_site)

Bases: django.contrib.admin.options.ModelAdmin

Defines the admin interface for Licences.

Currently set as default.

`media`

INDICES AND TABLES

- *Index*
- *Module Index*
- *Search Page*

MODULE INDEX

A

- animal, 18
- animal.admin, 22
- animal.forms, 20
- animal.models, 19
- animal.urls, 22
- animal.views, 20

D

- data, 13
- data.admin, 18
- data.forms, 17
- data.models, 13
- data.urls, 18
- data.views, 18

G

- groups, 24
- groups.admin, 25
- groups.models, 25
- groups.views, 25

T

- timed_mating, 22
- timed_mating.admin, 24
- timed_mating.forms, 23
- timed_mating.models, 23
- timed_mating.urls, 24
- timed_mating.views, 24

INDEX

A

add_measurement (in module data.views), 18
Animal (class in animal.models), 19
animal (data.models.Measurement attribute), 15
animal (module), 18
animal.admin (module), 22
animal.forms (module), 20
animal.models (module), 19
animal.urls (module), 22
animal.views (module), 20
animal_detail (in module animal.views), 21
AnimalAdmin (class in animal.admin), 22
AnimalForm (class in animal.forms), 20
AnimalForm.Media (class in animal.forms), 20
AnimalForm.Meta (class in animal.forms), 20
AnimalInline (class in animal.admin), 22
animals (data.models.Experiment attribute), 14
animals (data.models.Treatment attribute), 16
Assay (class in data.models), 13
assay (data.models.Measurement attribute), 15
Assay.DoesNotExist, 13
Assay.MultipleObjectsReturned, 13

B

Breeding (class in animal.models), 19
breeding (in module animal.views), 21
Breeding (timed_mating.models.PlugEvents attribute), 23
breeding_all (in module animal.views), 21
breeding_change (in module animal.views), 21
breeding_detail (in module animal.views), 21
breeding_plugevent (in module timed_mating.views), 24
breeding_pups (in module animal.views), 21
BreedingAdmin (class in animal.admin), 22
BreedingForm (class in animal.forms), 20
BreedingForm.Meta (class in animal.forms), 20
BreedingPlugForm (class in timed_mating.forms), 23
BreedingPlugForm.Meta (class in timed_mating.forms), 24

C

Cage (class in animal.models), 20

change_plugevents (in module timed_mating.urls), 24
contact (data.models.Environment attribute), 14
create_plugevents (in module timed_mating.urls), 24

D

data (module), 13
data.admin (module), 18
data.forms (module), 17
data.models (module), 13
data.urls (module), 18
data.views (module), 18
delete_plugevents (in module timed_mating.urls), 24
Diet (class in data.models), 14
diet (data.models.Treatment attribute), 16
Diet.DoesNotExist, 14
Diet.MultipleObjectsReturned, 14
diet_set (data.models.Vendor attribute), 17
donor (data.models.Transplantation attribute), 16

E

Environment (class in data.models), 14
environment (data.models.Treatment attribute), 16
Environment.DoesNotExist, 14
Environment.MultipleObjectsReturned, 14
environment_set (data.models.Researcher attribute), 15
Experiment (class in data.models), 14
experiment (data.models.Measurement attribute), 15
Experiment.DoesNotExist, 14
Experiment.MultipleObjectsReturned, 14
experiment_detail (in module data.views), 18
experiment_detail_all (in module data.views), 18
experiment_list (in module data.views), 18
experiment_set (data.models.Researcher attribute), 15
experiment_set (data.models.Study attribute), 16
ExperimentForm (class in data.forms), 17
ExperimentForm.Meta (class in data.forms), 17

G

get_absolute_url() (data.models.Experiment method), 14
get_absolute_url() (data.models.Researcher method), 15
get_absolute_url() (data.models.Study method), 16

- get_absolute_url() (data.models.Treatment method), 16
 - get_absolute_url() (timed_mating.models.PlugEvents method), 23
 - get_contact_title_display() (groups.models.Group method), 25
 - get_feeding_state_display() (data.models.Experiment method), 14
 - get_injection_display() (data.models.Experiment method), 14
 - get_mode_display() (data.models.Pharmaceutical method), 15
 - get_next_by_date() (data.models.Experiment method), 14
 - get_next_by_PlugDate() (timed_mating.models.PlugEvents method), 23
 - get_next_by_transplant_date() (data.models.Transplantation method), 16
 - get_ordering_display() (data.models.Vendor method), 17
 - get_previous_by_date() (data.models.Experiment method), 14
 - get_previous_by_PlugDate() (timed_mating.models.PlugEvents method), 23
 - get_previous_by_transplant_date() (data.models.Transplantation method), 16
 - Group (class in groups.models), 25
 - Group.DoesNotExist, 25
 - Group.MultipleObjectsReturned, 25
 - group_set (groups.models.License attribute), 25
 - GroupAdmin (class in groups.admin), 25
 - groups (module), 24
 - groups.admin (module), 25
 - groups.models (module), 25
 - groups.views (module), 25
- I**
- Implantation (class in data.models), 14
 - implantation (data.models.Treatment attribute), 16
 - Implantation.DoesNotExist, 14
 - Implantation.MultipleObjectsReturned, 14
 - implantation_set (data.models.Researcher attribute), 15
 - implantation_set (data.models.Vendor attribute), 17
- L**
- License (class in groups.models), 25
 - license (groups.models.Group attribute), 25
 - License.DoesNotExist, 25
 - License.MultipleObjectsReturned, 25
 - LicenseAdmin (class in groups.admin), 25
 - limited_object_detail (in module timed_mating.urls), 24
 - limited_object_list (in module timed_mating.urls), 24
- M**
- mark_deactivated() (animal.admin.BreedingAdmin method), 22
 - mark_sacrificed() (animal.admin.AnimalAdmin method), 22
 - Measurement (class in data.models), 15
 - Measurement.DoesNotExist, 15
 - Measurement.MultipleObjectsReturned, 15
 - measurement_set (data.models.Assay attribute), 14
 - measurement_set (data.models.Experiment attribute), 14
 - MeasurementForm (class in data.forms), 17
 - MeasurementForm.Meta (class in data.forms), 17
 - media (animal.admin.AnimalAdmin attribute), 22
 - media (animal.admin.BreedingAdmin attribute), 22
 - media (animal.admin.StrainAdmin attribute), 22
 - media (animal.forms.AnimalForm attribute), 20
 - media (animal.forms.BreedingForm attribute), 20
 - media (data.forms.ExperimentForm attribute), 17
 - media (data.forms.MeasurementForm attribute), 17
 - media (data.forms.StudyExperimentForm attribute), 17
 - media (data.forms.StudyForm attribute), 18
 - media (data.forms.TreatmentForm attribute), 18
 - media (groups.admin.GroupAdmin attribute), 25
 - media (groups.admin.LicenseAdmin attribute), 25
 - media (timed_mating.admin.PlugEventsAdmin attribute), 24
 - media (timed_mating.forms.BreedingPlugForm attribute), 24
 - model (animal.admin.AnimalInline attribute), 22
 - model (animal.forms.AnimalForm.Meta attribute), 20
 - model (animal.forms.BreedingForm.Meta attribute), 20
 - model (data.forms.ExperimentForm.Meta attribute), 17
 - model (data.forms.MeasurementForm.Meta attribute), 17
 - model (data.forms.StudyExperimentForm.Meta attribute), 17
 - model (data.forms.StudyForm.Meta attribute), 18
 - model (data.forms.TreatmentForm.Meta attribute), 18
 - model (timed_mating.forms.BreedingPlugForm.Meta attribute), 24
- P**
- Pharmaceutical (class in data.models), 15
 - pharmaceutical (data.models.Treatment attribute), 16
 - Pharmaceutical.DoesNotExist, 15
 - Pharmaceutical.MultipleObjectsReturned, 15
 - pharmaceutical_set (data.models.Vendor attribute), 17
 - PlugEvents (class in timed_mating.models), 23
 - PlugEvents.DoesNotExist, 23
 - PlugEvents.MultipleObjectsReturned, 23
 - PlugEventsAdmin (class in timed_mating.admin), 24
 - PlugFemale (timed_mating.models.PlugEvents attribute), 23
 - PlugMale (timed_mating.models.PlugEvents attribute), 23
- R**
- Researcher (class in data.models), 15

Researcher (timed_mating.models.PlugEvents attribute), 23
 Researcher.DoesNotExist, 15
 Researcher.MultipleObjectsReturned, 15
 researchers (data.models.Experiment attribute), 14
 researchers (data.models.Treatment attribute), 16

S

save() (timed_mating.models.PlugEvents method), 23
 Strain (class in animal.models), 19
 strain (data.models.Study attribute), 16
 strain_detail (in module animal.views), 21
 strain_detail_all (in module animal.views), 21
 strain_list (in module animal.views), 21
 StrainAdmin (class in animal.admin), 22
 Study (class in data.models), 15
 study (data.models.Experiment attribute), 14
 study (data.models.Treatment attribute), 16
 Study.DoesNotExist, 16
 Study.MultipleObjectsReturned, 16
 study_experiment (in module data.views), 18
 StudyExperimentForm (class in data.forms), 17
 StudyExperimentForm.Meta (class in data.forms), 17
 StudyForm (class in data.forms), 17
 StudyForm.Meta (class in data.forms), 18
 surgeon (data.models.Implantation attribute), 15
 surgeon (data.models.Transplantation attribute), 16

T

timed_mating (module), 22
 timed_mating.admin (module), 24
 timed_mating.forms (module), 23
 timed_mating.models (module), 23
 timed_mating.urls (module), 24
 timed_mating.views (module), 24
 Transplantation (class in data.models), 16
 transplantation (data.models.Treatment attribute), 16
 Transplantation.DoesNotExist, 16
 Transplantation.MultipleObjectsReturned, 16
 transplantation_set (data.models.Researcher attribute), 15
 Treatment (class in data.models), 16
 Treatment.DoesNotExist, 16
 Treatment.MultipleObjectsReturned, 16
 treatment_set (data.models.Diet attribute), 14
 treatment_set (data.models.Environment attribute), 14
 treatment_set (data.models.Implantation attribute), 15
 treatment_set (data.models.Pharmaceutical attribute), 15
 treatment_set (data.models.Researcher attribute), 15
 treatment_set (data.models.Study attribute), 16
 treatment_set (data.models.Transplantation attribute), 16
 TreatmentForm (class in data.forms), 18
 TreatmentForm.Meta (class in data.forms), 18

V

Vendor (class in data.models), 16
 vendor (data.models.Diet attribute), 14
 vendor (data.models.Implantation attribute), 15
 vendor (data.models.Pharmaceutical attribute), 15
 Vendor.DoesNotExist, 17
 Vendor.MultipleObjectsReturned, 17